
USER MANUAL **BENUTZERHANDBUCH**

MQ50 Oxygen Concentrator | Edition: V1.1

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MQ50 Oxygen Concentrator | Edition: V1.1

Table of Contents

SPECIAL NOTES	3
SAFETY NOTICE	4
Chapter 1: Product Introduction	7
● Working Principle	8
● Flow Diagram	8
● Product Characteristic.....	8
Chapter 2: Conditions of Use	9
Chapter 3: Intended Use and Contraindication	9
● Intended Purpose	9
● Intended User	9
● Contraindication.....	9
Chapter 4: Technical Specification	10
● Classification	10
● Main Performance	10
● Standard of Product Execution.....	11
● Product Specification	12
● Symbol Key	13
Chapter 5: Unpacking, Transportation and Storage	14
● Unpacking	14
● Transportation and Storage Environment.....	14
Chapter 6: Equipment Description.....	14
● Components	14
● Parts and Description.....	15
● Display and Description	16

Chapter 7: Operating Instructions.....	19
Chapter 8: Cleaning & Routine Maintenance.....	23
● Clean Device Enclosure.....	23
● Clean the Sponge Mesh Filter.....	23
● Clean Humidifier	24
● Disinfect Humidifier	25
● Filter Element Replacement.....	25
● Molecular Sieve Replacement.....	25
● Professional Maintenance Service	25
● Reuse of the Device.....	26
Chapter 9: Troubleshooting Guide.....	27
Chapter 10: EMC Information.....	29
Chapter 11: Maintenance	33
● Limited Warranty.....	33
● How to contact SYSMED	33
Chapter 12: Waste Disposal.....	34
Appendix A: Circuit Diagram	35
Appendix B: Packing List	35

SPECIAL NOTES

- Please read this manual carefully before using this product and save it for future reference. Do not use this product or approved accessories without fully reading and understanding these instructions. Injury or damage may result.
- If you cannot understand this manual (user manual), please contact your local dealer or after-sales service center. Our company provides consulting and training services. Make sure that the user understands this manual (user manual) and uses this product.
- This product is intended to supply oxygen, not for first aid or sustain life.
- For patients who are in urgent need of oxygen, critically ill patients, and continuous oxygen inhalation, when using this product, other oxygen supply devices (such as compressed oxygen cylinders or liquid oxygen systems) must be configured as a backup product. To prevent this product from malfunctioning, affecting the patient's oxygen intake.
- When the device malfunctions, the user should contact the product supplier or the manufacturer's after-sales service center in time, and must not disassemble and repair it.
- Severe carbon monoxide poisoning patients are banned.
- Nasal cannula is used exclusively by sole user.
- The product shall not be serviced or maintained while in use with a patient.
- If the product has been stored for 6 months by the supplier, the product may only be sold after passing the required inspection.

SAFETY NOTICE

Please read the following information carefully before operating the oxygen concentrator.

Warning

A warning indicates the possibility of injury to the user or the operator.

- **Smoking during oxygen therapy is dangerous and is likely to result in facial burns or death. Do not allow smoking or open flames within the same room as the oxygen concentrator or any oxygen-carrying accessories. If you smoke, you must always turn the oxygen concentrator off, remove the cannula and leave the room where either the cannula or mask or the oxygen concentrator is located. If unable to leave the room, you must wait 10 minutes after you have turned the oxygen concentrator off.**
- **Open flames during oxygen therapy are dangerous and are likely to result in fire or death. Do not allow open flames within 2 m of the oxygen concentrator or any oxygen-carrying accessories.**
- **There is a risk of fire associated with oxygen enrichment during oxygen therapy. Do not use the oxygen concentrator or accessories near sparks or open flames.**
- **To ensure receiving the therapeutic amount of oxygen delivery according to your medical condition, the oxygen concentrator must:**
 - 1) **be used with setting that have been individually determined or prescribed for you at your activity levels with your accessories;**
 - 2) **be used with the specific combination of parts and accessories that are in line with the specification of the concentrator or accessory manufacturer.**
- **The use of oxygen concentrator beyond the range of operating temperature, relative humidity or atmospheric pressure (altitude) specified in the conditions of use will affect the flow and oxygen concentration, and further affect the quality of the therapy.**
- **Pay attention to reduce fire risk when using oxygen therapy. Any material that is flammable or non-flammable in the air becomes extremely flammable and burns quickly at very high oxygen concentrations. For safety reasons, all kindling should be kept away from the oxygen concentrator, preferably not placed in the same room as the oxygen concentrator.**
- **Use only water-based lotions or salves that are oxygen-compatible prior to and during oxygen therapy. Never use petroleum-based or oil-based lotions or salves to avoid the risk of fire and burns.**
- **Do not lubricate fittings, connections, tubing, or other accessories of the oxygen concentrator to avoid the risk of fire and burns. Only use approved lubricants.**
- **Periodically evaluate the therapeutic effect of the device, contact your doctor if you feel not good or uncomfortable.**
- **When using the oxygen concentrator, please follow your doctor's instructions and the**

requirements of the user manual. If the patient or caregiver notices that oxygen supply is insufficient, they should contact the product supplier or healthcare provider immediately. Do not adjust oxygen flow unless directed by a healthcare provider.


- Oxygen makes it easier for a fire to start and spread. Do not leave the nasal cannula or mask on bed coverings or chair cushions, if the oxygen concentrator is turned on, but not in use; the oxygen will make the materials more flammable. Turn the oxygen concentrator off when not in use to prevent oxygen enrichment.
- Do not use it when showering. If the patient needs to use the device continuously, the device must be placed at least 3 m away from the bathroom.
- Don't touch the oxygen concentrator when your body is wet. Do not use or store the device near water or other fluids that easily conduct electricity.
- Do not contact the oxygen concentrator with water or other easily conductive liquids. If it falls into water, unplug it immediately.
- Small parts that come off the equipment, such as ingestion or inhalation, can cause asphyxia.
- Use only spare parts recommended by the manufacturer to ensure proper function and to avoid the risk of fire and burns.
- Do not pull or stretch the power cord to move the device.
- Do not drop any foreign matter into any opening.
- Modification of this product is strictly prohibited.
- Geriatric, paediatric or any other patient unable to communicate discomfort can require additional monitoring and or a distributed alarm system to convey the information about the discomfort and or the medical urgency to the responsible care giver to avoid harm.
- If you feel discomfort or are experiencing a medical emergency while undergoing oxygen therapy, seek medical assistance immediately to avoid harm.
- Do not connect this product in parallel or in series with other brands of oxygen concentrators or oxygen therapy equipment.
- The oxygen concentrator cannot be used in environments, such as electrocautery, electrosurgery, defibrillation, X ray (γ ray), Infrared radiation, transient electromagnetic fields, including magnetic resonance (MRI) and radio interference.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be more than 30cm away from any part of this oxygen concentrator. Otherwise, the performance may be reduced.

Cautions

- When using it, the oxygen concentrator will produce regular exhaust sound, which is a normal phenomenon.
- In order to achieve the best performance, it is recommended that the oxygen concentrator should be turned on for more than 30 minutes each time. Frequent operation in a short time will reduce the service life of the oxygen concentrator.
- The air inlet and exhaust ports of the oxygen concentrator should be located in a well-ventilated place. The air inlet should be kept away from plush, hair or other similar objects, and it is strictly forbidden to block the air inlet of the oxygen concentrator; It is strictly forbidden to place the oxygen concentrator on a soft surface, such as a sofa or bed, where the exhaust port can be blocked.
- The oxygen concentrator should be at least 40 cm away from walls, tapestries, furniture etc.
- Each medical oxygen concentrator is equipped with 2 sets of nasal cannula before delivery. Patients can also choose other anti-extrusion registered nasal cannula according to their needs, but no longer than 4m.
- The use of unspecified humidifier, nasal cannula or other accessories may reduce the use effect of the oxygen concentrator.
- Be careful not to place the oxygen concentrator in an environment that is prone to collision or tripping over and entanglement by power cords and nasal cannula. Power cords and nasal cannula pose a risk of strangulation caused by entanglement in infants or children.
- Unexpected movements of children and pets may affect the normal operation of the equipment system.
- The oxygen concentrator, its parts and accessories are specified for use at specific flows.
- The oxygen concentrator and its components or accessories are compatible.
- The PATIENT is an intended OPERATOR.
- This product may cause radio interference or disrupt the operation of nearby equipment. Mitigation measures may be necessary, such as reorienting and positioning the equipment or shielding the site.
- Most appliances are susceptible to radio frequency interference. Therefore, using portable communication equipment near the oxygen concentrator may cause interference to the machine.
- The oxygen concentrator should not be stacked or put too near other devices. This may result in improper operation. If it has to be like that, pay attention to observe and inspect the status of the oxygen concentrator to ensure it works normally.

Chapter 1: Product Introduction


Welcome to use our SYSMED medical oxygen concentrator!

SysMed (China) Co., Ltd. was invested by Chinese Academy of Science, Shenyang Institute of Automation and Dalian Institute of Chemical Physics. We have been committed to professional oxygen concentrator quality and professional oxygen therapy service to requite your love for SYSMED.

● Working Principle

The oxygen concentrator uses the principle of Pressure Swing Adsorption (PSA) of molecular sieve to obtain oxygen by adsorbing nitrogen and other gas components in the air.

● Flow Diagram

● Product Characteristic

SYSMED medical oxygen concentrator guarantees long-term oxygen supply with stable oxygen concentration and reliable operation with 6 professional measures.

- Use oil-free compressor with large displacement, sufficient raw gas; use U.S. wear-resistant seal-ring, less wear and tear, large gas output.
- Using molecular sieve automatic filling technology, pre-loading force is large, not easy to crush; The molecular sieve tank adopts aluminum profile, automatic pressing, high fatigue strength, and maximizes the service life.
- Good heat dissipation, 24 hours continuous running, stable oxygen concentration.
- Using independently developed oxygen concentration monitoring technology (invention patent No. ZL200810011572.4), first in China.
- Using unchangeable cumulative time function to indicate the device's status, provide basis for regular testing (inspection) and maintenance service, and provide truly reliable basis for fulfilling quality warranty promise.
- Multi-variable control technology of time and pressure can be used to adjust system parameters through maintenance service, so that the device performance is always in the best.

Chapter 2: Conditions of Use

- Electrical Requirements: 220/ 230V~, 50Hz
- Operating Temperature: 10°C ~ 37°C
- Operating Relative Humidity: 20% ~ 65%
- Operating Atmospheric Pressure: 80kPa ~ 101kPa
- The surrounding environment is dry and ventilated, free from pollutants or fumes, and has no strong electromagnetic interference.

Caution:

- When the AC power supply voltage exceeds -15% to +10% of the rated voltage, the voltage is too high, which may result in damage to the equipment. If the voltage is too low, the equipment may not start. If the grid is unstable, please install a voltage regulator before use.
- It is recommended for use in terrestrial environments below 1828m above sea level. In the environment above 1828m above sea level, it may lead to a decrease in the adsorption efficiency of the oxygen concentrator and a decrease in oxygen concentration.
- The main control board and electronic components may fail if the ambient humidity is too high.
- Below 10°C, compressor start-up may be difficult. Above 37°C, the compressor may overheat and shorten the service life of the compressor.

Chapter 3: Intended Use and Contraindication

● **Intended Purpose**

Oxygen concentrator is intended for the administration of supplemental oxygen. The device is not intended for life support nor does it provide any patient monitoring capabilities.

● **Intended User**

This product is suitable for people with insufficient oxygen supply, or people with weak body and poor body immunity. When using, follow the doctor's instructions and select the applicable parameters. Infants should not use this product.

● **Contraindication**

The product is prohibited for those who have oxygen poisoning and oxygen allergy.

Chapter 4: Technical Specification

● Classification

- Type of protection against electric shock: class II
- Degree of protection against electric shock: Type BF Applied Part

NOTE: The nasal cannula is the applied part of the oxygen concentrator.


- IP classification: IP21
- The product must NOT be used in the presence of flammable anesthetic mixtures with air or with oxygen or nitrous oxide
- Mode of operation: Continuous operation

● Main Performance

- Within 30 min of starting the product, the oxygen output and oxygen concentration can reach the rated requirements.
- Oxygen outlet output pressure: $0.05\text{MPa} \pm 0.005\text{MPa}$ (in normal condition)
 $<0.15\text{MPa}$ (in single fault condition)


NOTE: The measurement uncertainty of oxygen outlet output pressure is 0.003MPa.

- The pressure release range of oxygen concentrator is 250 kPa-270 kPa.
- The expected service life of the product: 5 years (under normal working environment and normal maintenance conditions).
- Production date: See device label.
- Under STPD (standard temperature and atmospheric pressure, dry) conditions, the rated range of oxygen output flow and corresponding oxygen concentration is shown in the figure below:

Relationship of Outlet Concentration and Flow(STPD)

- Under the conditions of different altitudes (atmospheric pressure), the rated range of oxygen output flow and corresponding oxygen concentration is shown in the following figure:

● Standard of Product Execution

Standard No.	Standard Title
IEC 60601-1: 2020	Medical electrical equipment-Part 1: General requirements for basic safety and essential performance
IEC 60601-1-2: 2020	Medical electrical equipment-Part 1-2: General requirements for basic safety and essential performance-Collateral Standard: Electromagnetic disturbances -Requirements and tests
IEC 60601-1-8: 2020	Medical electrical equipment-Part 1-8: General requirements for basic safety and essential performance-Collateral Standard: General requirements, tests and guidance for alarm systems in medical electrical equipment and medical electrical systems
IEC 60601-1-11: 2020	Medical electrical equipment-Part 1-11: General requirements for basic safety and essential performance-Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment
ISO 80601-2-69: 2020	Medical electrical equipment- Part 2-69: Particular requirements for the basic safety and essential performance of oxygen concentrator equipment

● Product Specification

Model	MQ50	
Rated Power (VA)	320	
Standby Power (VA)	0	
Rated Voltage	220/ 230V~, 50Hz	
Flow rate(l/min)	1~2	3~5
Concentration (V/V)	>90%	>87%
Sound pressure level dB(A) (Adjustable flow range)	<40	
Sound power level dB(A) (Adjustable flow range)	<48	
Net weight (kg)	15.5	
Dimension(mm)	390×230×600	
Function	System high temperature alarm; Low oxygen concentration alarm; Power failure alarm; Pressure abnormal alarm; Low flow alarm; NTC sensor error alarm; Timer.	

NOTE:

- The flow rate specifications listed in the table are expressed at STPD (standard temperature and pressure, dry);
- The measurement uncertainty of sound pressure level and sound power level is 2dB(A);
- The measurement uncertainty of oxygen concentration is included in the values listed in the table.

● Symbol Key

Graphic	Title	Graphic	Title
	"ON" (power)		"OFF" (power)
	Follow Instruction for Use		General warning sign
	CLASS II equipment		Type BF Applied Part
	CE certification mark		Alternating current
	No smoking		No open flame: Fire, open ignition source and smoking prohibited
	Operating instructions		No sitting
	Output; exit		Do not reuse
	Warning, electricity		Authorized Representative in the European Community
	Manufacturer		Catalogue number
	Batch code		Serial Number
	Use by date		Date of Manufacture
	Stacking Limit By Number		This Way Up
	Fragile, handle with care		Keep dry
	Temperature limit		Atmospheric pressure limitation
	Humidity limitation		Do not use if package is damaged
	Unique device identifier		Medical device
IP21	Against ingress of solid foreign objects: ≥ 12.5 mm diameter Against ingress of water with harmful effects: vertically dripping		WEEE-Subject to waste electrical and electronic equipment regulations, i.e. not for general waste

Chapter 5: Unpacking, Transportation and Storage

● Unpacking

- Check whether the carton is damaged or not placed vertically. If damage is evident, inform the carrier and the product supplier.
- Take out the protective foam block and accessories from the carton.
- Carefully take out the concentrator and other related parts from the carton. Inspect/examine the exterior of the oxygen concentrator for nicks, dents, scratches or other damage. Check all components whether as listed. Please contact the product supplier or manufacturer's after-sales service center if there are short of components or quality issues.

Caution:

- Please keep the carton and packing materials for storage and transport.

● Transportation and Storage Environment


- This product must be placed upright and vertical during transportation. It is strictly prohibited to place it horizontally or upside down.
- Transport & Storage Temperature: $-30^{\circ}\text{C} \sim 70^{\circ}\text{C}$.
- Transport & Storage Relative Humidity: 15%~95%, no condensation.
- Transport & Storage Atmospheric Pressure: 500hPa~1060hPa.

Chapter 6: Equipment Description


● Components

The oxygen concentrator is composed of compressor, filter, molecular sieve adsorption separation device, control device, flow meter, humidifier, nasal cannula and so on.

● Parts and Description

Front View of Device

Back View of Device

No.	Parts	Description
1	Flowmeter	The position of the float indicates the output oxygen flow rate
2	Oxygen Outlet	Oxygen output
3	Serial Port	The serial port is used only for internal debugging of the manufacturer. No operation is allowed on the serial port by non-internal personnel
4	Flowmeter Adjustment Dial	Adjust the output oxygen flow
5	Humidifier	Connection thread specification: M14 Volume: 350ml Humidifier is used to moisten oxygen and prevent dry oxygen from irritating the throat and nasal mucosa. As well as to prevent phlegm from drying up and not being spit out. It is composed of a humidifier bottle, a humidifier cap and a connector.
6	LCD	Displays product operation parameters and information.
7	Power LED	Indicates the device is running. It is green when it lights up.
8	Alarm LED	Indicates the alarm state.
9	Power Switch	Power on/ Power off


10	Model Label	Product model identification
11	Power Cord	Connect power supply.
12	Circuit Breaker	Over-current protection
13	Air Intake Window	Open to replace the sponge mesh filter.
14	Humidifier Band	Fix the humidifier to the device.
15	Filter Access Door	Open to replace the inlet filter.
16	Device Label	Product information identification.

● Display and Description

- Startup display interface (for about 5 seconds)


- Working display interface

NOTE: The high oxygen icon "hoch" flashes within 5 minutes after startup.

The device has an audible alarm and LED indicators, as shown above LCD Screen instructions:

LCD Zone	Indicators	Possible Cause	Status Prompt	Test Method	Alarm delay time	Priority
①		The audio alarm	----	Visual inspection	----	----
②	Hohe Temp	Compressor temperature is too high	Red LED flashes twice per second and the Audible Alarms is sounding hurriedly. The device is not operating.	Can only be operated by a technical engineer authorized by SYSMED	<1 min	High
	Druck abnormal	Working pressure abnormal	Red LED flashes twice per second and the Audible Alarms is sounding hurriedly. The device is not operating.	Can only be operated by a technical engineer authorized by SYSMED	<1 min	High
	NTC Sensor Fehler	Temperature sensor failure	Yellow LED flashes once per 2 seconds and the Audible Alarms is sounding intermittently.	Can only be operated by a technical engineer authorized by SYSMED	<1 min	Medium
	Niedriger Flow	Output flow is too low	Yellow LED flashes once per 2 seconds and the Audible Alarms is sounding intermittently.	Can only be operated by a technical engineer authorized by SYSMED	<1 min	Medium
③		High Oxygen Concentration	----	Visual inspection	----	----

		Low Oxygen Concentration Alarm: Oxygen concentration below 82%	Red LED flashes twice per second and the Audible Alarms is sounding hurriedly.	Adjust the oxygen flow-meter to the maximum level until alarm	<1min	High
④	"00:00"	Current working time	----	Visual inspection	----	----
⑤	"00000"	Elapsed time (Hours)	----	Visual inspection	----	----
		This symbol displays when alarm occurs	----	Visual inspection	----	----
	----	Power outage alarm	Audible Alarm is sounding continuously, and alarm duration >60 seconds.	Unplug the power plug while in operation	----	----

 Caution:


- All of the above alarms are technical alarms.
- When an alarm occurs, the operator can identify the failure type from the LCD and contact homecare provider immediately.
- Operation position: Should be in any position that is vertical or horizontal at a 30° angle to the plane of the LCD display.
- The recommended test interval for the alarm system is: 18 months.
- Medium priority alarms will be accompanied by an auditory alarm signal of "Du.. Du.. Du". High priority alarms will be accompanied by an auditory alarm signal of "DuDuDu...DuDu, DuDuDu...DuDu".
- When the equipment is switched on, the alarm system will be self-checked. If the alarm system is in normal operation, a "Du" sound will be generated to prompt.
- The alarm settings before the power loss will not be affected, regardless of the duration of the power loss.
- All variables for display, control, and alarm are calculated using arithmetic average filtering method.
- The alarm sound pressure level is 50-70 dB(A). If the alarm volume is not within this range, please contact the product supplier or the manufacturer's after-sales service center.

Chapter 7: Operating Instructions

1. The oxygen concentrator should be placed on flat ground in a clean, ventilated and direct sunlight-free environment, and keep the oxygen concentrator at least 40cm away from walls and other objects.

Warning:

- Do Not put sundries, water and oil containers and other covers on the oxygen concentrator.
 - Do Not place any sundries at the bottom of the oxygen concentrator, so as not to block the exhaust outlet at the bottom, causing overheating and affecting normal operation.
2. After reading this entire manual, plug the power cord into an electrical outlet.
3. Unscrew the humidifier cap in a counterclockwise direction, remove the transparent humidifier bottle, fill with a proper amount of distilled water, and then tighten the brown cap clockwise.


Caution:


- Please use distilled water, purified water or cold boiled water for the humidifier, and replace it once a day.
 - The level of water should be between the min and the max markings. Filling too much may overflow.
4. Install the humidifier with water in the hook and loop fastener on the top of the device. Use the hook and loop fastener to tighten the bottle and secure the bottle.

5. Please connect the PVC hose to the oxygen concentrator outlet and humidifier inlet.

6. Please connect the fire damper to the outlet end of the humidifier.


⚠ Caution:

- One end of the fire damper is an air inlet and the other end is an air outlet. Pay attention to the gas flow sign on the fire damper, and do not change the installation direction of the fire damper. Do not remove the connection tube from the fire damper during installation.

7. Please connect the intake end of the nasal cannula to the air outlet of the fire damper.

8. The air path connection, see the Removable Parts Connection Diagram below.

Removable Parts Connection Diagram

⚠ Warning:

- This product must be connected to the fire damper according to the above instructions, otherwise SYSMED will not be responsible for any consequences!

9. Turn on the power switch.

⚠ Caution:

- The power cord should be kept away from heat sources or high-temperature surfaces.
- When the storage temperature is the lowest storage temperature or the highest storage temperature, the time required for the oxygen concentrator to be ready for its intended use when the ambient temperature is 20 °C is 4 hours.

10. Block the outlet of the humidifier with your finger for about 20 seconds. Observe whether the float in the flowmeter drops, listen if there is a " Chu, Chu, Chu " sound from the humidifier. If so, it indicates that the humidifier is tight. If not, it indicates the humidifier leaks. Take off the humidifier, remove and tighten the cap again. If there is still no sound after adjustment, please contact the product supplier or the manufacturer's after-sales service center.

11. Adjust the oxygen output flow with eyes, scale line and the center of the float on the same line.

⚠Caution:

- The flowmeter adjustment dial increases the flow counterclockwise and decreases the flow clockwise.

⚠Warning:

- Please do not adjust the position of the float of the flowmeter beyond the red line, otherwise it will reduce the efficiency of the oxygen concentrator in long-term use.
- Home users should strictly follow the doctor's instructions when doing oxygen therapy. Oxygen therapy should be carried out according to the oxygen flow rate set by the doctors and the time of daily oxygen flow inhalation. It is strictly forbidden to adjust oxygen flow and the length of time.
- Do not turn on the oxygen concentrator when the flow-meter is off.

12. Insert the nasal plug on the nasal cannula into the nostrils of the patient for oxygen inhalation.

⚠Caution:

- When inserting and removing the nasal cannula, apply moderate force to achieve the sealing effect.
- Reach in front of the nasal cannula outlet, and if you can't feel the airflow, check the nasal cannula connection for air leaks, or you can place the nasal cannula outlet under water to observe for bubbles.
- The outer diameter of the outlet of the oxygen concentrator and the outlet of the humidifier connected to the nasal cannula is $\Phi 6\text{mm}$.
- Users should purchase a nasal cannula with a medical device registration certificate and a valid period. For specific use requirements, please refer to the instructions for the purchased nasal cannula.
- The position of the nasal cannula in the nostrils and the orientation of its air outlet determine the amount of oxygen output to the user's respiratory system.
- The nasal cannula should be placed in a good way to prevent the danger of getting entangled around the head or neck.
- The supplied nasal cannula is a single-use product, which is easy to breed bacteria and cause cross-infection after long-term use.

13. When the oxygen therapy is finished, turn off the power switch.

⚠Caution:

- Unplug the power cord when it is not used continuously.
- Please don't move the oxygen concentrator when using.

Chapter 8: Cleaning & Routine Maintenance

⚠Caution:

- In normal use, nasal cannula and humidifier may be contaminated by bodily fluids or respiratory gases.
- The cleaning and maintenance of nasal cannula should refer to its user manual.
- If daily use of the oxygen concentrator is not required per medical requirements or doctor's instructions, operate the device for a minimum of 1 hour weekly to maintain performance and prolong its service life.

⚠Warning:

- Unplug the power cord before cleaning and maintenance to prevent electric shock.
- Do NOT immerse the product in any liquid.
- Do NOT use alcohol, acetone, chlorides or abrasive cleaners.
- Ensure no residue of the detergent is left after cleaning them with a neutral detergent. Residue of the detergent in the airway may cause damage to the patient's respiratory system.
- Check whether the product and its components are damaged after cleaning. If damaged, please contact the product supplier promptly for replacement.

● Clean Device Enclosure


1. Use a damp cloth with warm water and neutral detergent to gently wipe the device enclosure.
2. Completely remove the residue of the detergent with a damp cloth.
3. Allow the concentrator to air dry, or use a dry cloth, before operating the concentrator.

⚠Caution:

- Ensure the device is completely dry before operating it.
- It is recommended to clean it once or twice per month.

● Clean the Sponge Mesh Filter

1. Open the air intake window, remove the sponge mesh filter and dust it off.

2. Gently wash the sponge mesh filter with warm water and neutral detergent, then rinse it thoroughly.
3. Dry it naturally thoroughly in a clean and cool place.
4. Reinstall the sponge mesh filter. Or a spare air intake filter can be used alternately.

⚠ Caution:

- Do NOT vigorously rub the filter during the cleaning process.
- Do NOT expose the sponge mesh filter to strong sunlight for a long time.
- Before operating the concentrator, make sure that the sponge mesh filter has been put in the air intake window.
- Replace the sponge mesh filter if any damage is found.
- It is recommended to clean weekly to prevent the growth of bacteria. If the environment contains a large amount of dust or smoke, more frequent cleaning and replacement may be necessary. Delayed replacement may affect your health.

● Clean Humidifier

1. Open the cap of the humidifier, clean the humidifier thoroughly with warm water and neutral detergent, and then rinse it.
2. Wipe it with a clean, dry cloth, and air dry thoroughly.
3. Wrapped the humidifier in a clean plastic bag after cleaning when not in use.

⚠ Caution:

- Care is required when cleaning the humidifier to avoid bumps.
- Do NOT use hard or abrasive cleaning tools to brush the humidifier to avoid scratching the surface.
- It is recommended to clean it once a day to prevent the growth of mold and bacteria.
- The life expectancy of the humidifier is 5 years. Please replace the humidifier if there is a decrease in flow due to mineral deposition.

● Disinfect Humidifier

If you clean your humidifier according to the above cleaning requirements during use, it is not necessary to disinfect your humidifier. When the humidifier is contaminated or used by different users, it should be disinfected according to the disinfection method described below.

1. Clean the humidifier thoroughly according to the above method.
2. Soak the cleaned humidifier in the disinfectant solution for 20 to 30 minutes and rinse it with warm or hot water.
3. Wipe it with a clean, dry cloth and dry naturally.

⚠Caution:

- Generally, the disinfectant that can be used in the home is 1:1 white vinegar and distilled water, or 1:10 bleach and distilled water.

● Filter Element Replacement

When the oxygen concentrator is used in the indoor environment, it is recommended to replace the inlet filter after running for every 4,000 to 8,000 hours, depending on environmental factors (e.g., air pollution, humidity). Please contact the product supplier or the after-sales service center of the manufacturer.

1. Open the filter access door and pull out the inlet filter.

2. Reinstall a new inlet filter. Then close the filter access door.

● Molecular Sieve Replacement

Under rated operating environment and normal working condition, it is recommended that the replacement cycle of molecular sieves is 15,000 hours.

● Professional Maintenance Service

SYSMED medical oxygen concentrators are designed to minimize routine maintenance. Only professionals in the manufacturer's after-sales service center or people who are very familiar

with the operation process, such as authorized personnel or factory-trained staff, can carry out repairs or commissioning.

To order any accessories, please contact the product supplier.

If you need assistance when installing, using, or accessories, please contact your clinician or product supplier.

 Caution:

- The product and accessories should be inspected regularly or before use, and if damaged, please contact the product supplier to replace them in time.
- In addition to the above-mentioned routine maintenance work, other service work can only be carried out by technical engineers authorized by SYSMED. If the product is damaged by unauthorized repairs, the warranty period of the product is void and additional repair costs are required.
- In order to ensure that the product can be used for a long time, the user must comply with the safety instructions and cleaning and maintenance requirements of the user manual.

● Reuse of the Device

The MQ50 concentrator is suitable for reuse with other patients.

The concentrator must be cleaned and serviced before reuse. This work may only be carried out by the manufacturer or a third party authorized by the manufacturer.

The following points must be observed when reusing the concentrator:

- Dispose of all non-reusable accessories such as the oxygen tube, nasal cannula, and humidifier.

 **Caution: The concentrator must be disconnected from the power supply for the following steps!**

- Remove dust deposits from inside the device housing.
- Clean/disinfect the device housing inside and out. Use a suitable disinfectant such as Dismozon Plus or Microbac Forte. Follow the manufacturer's instructions for the disinfectant.
- Replace the sponge mesh filter and the inlet filter as described in Chapter 8.
- Replace all damaged and defective parts.
- Check the oxygen concentration.

Chapter 9: Troubleshooting Guide

Problem	Why it Happened	What to Do
The device is switched on but not working.	Internal part failure	Contact your product supplier or manufacturer's after-sales service center.
The device is switched on but not working, or only working intermittently.	The power cord is damaged.	Check whether the power cord is damaged.
		Contact your product supplier or manufacturer's after-sales service center.
	The power cord plug is not properly inserted into the electrical outlet.	Make sure the device is properly plugged into the electrical outlet.
	The oxygen concentrator is not receiving power from the electrical outlet.	Check your household outlet fuse or circuit.
		Contact your product supplier or manufacturer's after-sales service center.
Oxygen is not flowing or the oxygen flow is limited.	The nasal cannula is kinked or blocked, stopping the delivery of oxygen.	Unkink the nasal cannula or replace if necessary.
	The humidifier is not properly connected to the device.	Reinstall the humidifier.
The float in the flow meter does not move up or down when adjusting the flow control dial.	The flow control dial is not open.	Turn the flow control dial slowly and carefully.
	The flow control dial is faulty.	Contact your product supplier or manufacturer's after-sales service center.
There is water in the cannula	There has been a change in temperature; OR The device is too close to a wall, curtains or furniture.	Dry the inside of the humidifier cap.
		Do not use hot water in the humidifier.
		Do not overfill the humidifier.
		Keep the oxygen concentrator and cannula in the same room at the same temperature.

 Caution:

- For any product-related inquiries including device malfunctions, operational issues, or treatment application concerns, please immediately contact: your healthcare provider, the authorized distributor, the manufacturer's customer service center.
- We also welcome your feedback regarding user experience and product improvement suggestions through our official service channels.

Chapter 10: EMC Information

EMC refers to the ability of a product to suppress the electromagnetic interference (EMI) to other devices, while not causing similar electromagnetic radiation interference to other devices. The product may cause electromagnetic interference to other devices through the air or by connecting to a power cord.

To ensure the electromagnetic compatibility of the product, installation, commissioning, and use should be carried out according to the user manual.

Portable and mobile radio frequency (RF) communication devices may affect the electromagnetic compatibility of the product. If this occurs, please contact SYSMED's after-sales service center for solutions.

Solutions to common EMC problems:

- Strictly follow the operation procedures specified in the product user manual to ensure the basic working conditions of the product are not affected by electromagnetic interference.
- Keeping other devices away from this product can reduce electromagnetic interference.
- Adjusting the relative position or installation angle between this product and other devices can reduce electromagnetic interference.
- Changing the wiring location or signal cables location of other devices can reduce electromagnetic interference.
- Changing the power path of other devices can reduce electromagnetic interference.

This product is classified as Group I Class B equipment by CISPR 11 Standard, non-life support device.

Guidelines and Manufacturing's statements – Electromagnetic emission: the product is intended to be used in the following specific electromagnetic environments, and the user should ensure that the product is used in such an environment.

Emissions Test	Compliance	Electromagnetic environment - Guidance
RF emissions CISPR 11	Group I	The product uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The product is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/Flicker emissions IEC 61000-3-3	Complies	

Guideline and manufacturer's Statement – electromagnetic immunity:

The product is intended for use in the electromagnetic environments specified below, and the user of this product should make sure it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic environment guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±8kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical Fast Transient/Burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input-output lines	±2 kV for power supply lines Not applicable	Mains power quality should be that of a typical home or hospital environment
Surge IEC 61000-4-5	±0.5 kV, ±1 kV differential mode ±0.5 kV, ±1 kV, ±2 kV common mode	±0.5 kV, ±1 kV differential mode Not applicable	Mains power quality should be that of a typical home or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT; 1 cycle at 0° 70% UT; 25 cycles at 0° 0% UT; 250 cycles	0% UT; 1 cycle at 0° 70% UT; 25 cycles at 0° 0% UT; 250 cycles	Mains power quality should be that of a typical home or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical hospital or home environment.
NOTE: UT is the A.C. mains voltage prior to application of the test level			

Guidance and Manufacturer's declaration – Electromagnetic Immunity:

This product is intended for use in the electromagnetic environment specified below. The user of this product should make sure it is used in such an environment.

Immunity Test	IEC 60601 Test level	Compliance level	Electromagnetic environment-Guide
Conducted RF IEC 61000-4-6	3 V 150kHz-80MHz 6Vrms Amateur Radio & ISM Bands between 150 kHz and 80 MHz	3V 150kHz-80MHz 6Vrms Amateur Radio & ISM Bands between 150 kHz and 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d=1.2 \sqrt{P}$ $d=1.2 \sqrt{P} \quad 80\text{MHz}-800\text{MHz}$ $d=2.3 \sqrt{P} \quad 800\text{MHz}-2.5\text{GHz}$ P – where P is the maximum output power rating of the transmitter in Watts (W) according to the transmitter manufacturer. d-where d is recommended separation distance in meters (m) The field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol. 
Radiated RF IEC 61000-4-3	10 V/m 80MHz to 2.7GHz	10 V/m 80MHz to 2.7GHz	
NOTE 1: At 80MHz and 800MHz, the higher frequency range applies.			
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a: Field strengths from fixed transmitters, such as base station for radio (cellular/cordless) telephones and			

land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicated theoretically with accuracy. To access the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If an abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

b: Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

Recommended isolation distance between portable and mobile RF communication equipment and equipment or system: the product is intended to be used in a controlled electromagnetic environment of radiation radio frequency disturbance, the purchaser or user of the product can prevent electromagnetic interference by maintaining a minimum distance between the portable and mobile radio frequency communication device and the product, as recommended below, depending on the maximum output power of the communication device.

Rated maximum output power of transmitter (W)	Separation distance according to different frequencies of transmitter/m		
	150kHz~80MHz	80MHz-800MHz	800MHz-2.5GHz
	$d=1.2\sqrt{P}$	$d=1.2\sqrt{P}$	$d=2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in Watts (W) according to the transmitter manufacturer.

NOTE 1: At 80MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

The product has passed the electromagnetic compatibility standard test, in the electromagnetic interference environment, the basic performance of the product:

- Under electromagnetic interference, adjust flow-meter to rated flow, the product can create oxygen. When the voltage is reduced, the operation of the product is intermittent, when flow-meter at 50% of rated flow, oxygen purity meets the rated requirements. When voltage sag interruption, the product automatically restored to oxygen produce statement.
- Default state: The default statement of the product is factory setting.

Chapter 11: Maintenance

● Limited Warranty

Sysmed (China) Co., Ltd. warrants that the system shall be free from defects of workmanship and materials and will perform in accordance with the user manual for a period of three (3) years from the date of sale by Sysmed (China) Co., Ltd. to the dealer.

Sysmed (China) Co., Ltd. disclaims all liability for economic loss, loss of profits, overhead, or consequential damages which may be claimed to arise from any sale or use of this product. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

This warranty is given in lieu of all other express or implied warranties, including the implied warranties of merchantability and fitness for a particular purpose. In addition, in no event shall SYSMED be liable for lost profits, loss of good will, or incidental or consequential damages, even if SYSMED has been advised of the possibility of the same. Some states or provinces do not allow the exclusion or limitation of implied warranties or the disclaimer of incidental and consequential damages. Accordingly, the laws of your state or province may give you additional protections.

To exercise your rights under this warranty, contact your local authorized Sysmed (China) Co., Ltd. dealer or contact Sysmed (China) Co., Ltd.

● How to contact SYSMED

To have your device serviced, contact your product supplier. If you need to contact SYSMED directly, call the SYSMED After-Sales Service Center at: +8602423970166. You can also use the following address:

Sysmed (China) Co., Ltd.

No.299, Baita Third Street, Hunnan District, 110169 Shenyang, Liaoning,
PEOPLE'S REPUBLIC OF CHINA.

Visit the MQ web site at: <http://www.sysmed.cn/>

Chapter 12: Waste Disposal


Dispose of the device in accordance with local regulations.

WEEE/RoHS Recycling Directive

When the product is at the end of its life and the user intends to discard the product, it must be disposed of separately from other production waste. For electronic and electrical waste (e.g., circuit boards) and molecular sieve, do NOT dispose of it inappropriately. As well as discarded nasal cannula produced during the use of the product. Please contact your local agency or waste disposal service center.

Other waste such as filter, shells and other waste of the product have no impact on environmental protection and pollution, if they are disposed of in accordance with the relevant regulations of the local government.

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Corp.GmbH(Europe)
VibraCare GmbH